



2ml 9-425 Autosampler Vials & Closures

- 9-425 Vials are uniformly flat bottom for security with inserts.
- Pre-assembled caps and septa are convenient and minimize contamination from handling.



9-425 Vial



P/N	Size	Packing	Description
2ML-9-V1001	12x32	100PCS/PK	2ml, Clear Glass, Flat Base, 9-425 Screw Thread Vial
2ML-9-V1002	12x32	100PCS/PK	2ml, Clear Glass, Flat Base, 9-425 Screw Thread Vial, With label
2ML-9-V1003	12x32	100PCS/PK	2mL, Amber Glass, Flat Base, 9-425 Screw Thread Vial
2ML-9-V1004	12x32	100PCS/PK	2mL, Amber Glass, Flat Base, 9-425 Screw Thread Vial, With label
2ML-9-SV1001	12x32	100PCS/PK	0.3ml, 9-425 PP Screw Thread, Micro-Vial, Transparent
2ML-9-SV1002	12x32	100PCS/PK	0.3ml, 9-425 PP Screw Thread, Micro-Vial, Amber

05

9-425 Cap



Smooth type



Ribbed type

P/N	Size	Packing	Description
9-SP2001-1	9-425	100PCS/PK	Blue, 9-425 Open Top Smooth Screw Cap
9-SP2001-2	9-425	100PCS/PK	Blue, 9-425 Open Top Ribbed Screw Cap
9-SP2002	9-425	100PCS/PK	Black, 9-425 Open Top Screw Cap





9-425 PTFE/Silicone Septum

- PTFE/Silicone septa are most popular for HPLC applications, and Pre-slit septa are easier to pierce with needles.
- Mercedes shape slit septa has been put into use.



Pre-slit, Mercedes shape

P/N	Size	Packing	Description
9-SP1001	9-425	100PCS/PK	9mm Red PTFE/White Silicone 1mm(.040")thick
9-SP1002	9-425	100PCS/PK	9mm White PTFE/Red Silicone 1mm (.040 in) thick
9-SP1003	9-425	100PCS/PK	9mm Blue PTFE/White Silicone 1mm (.040 in) thick
9-SP1004	9-425	100PCS/PK	9mm Red PTFE/White Silicone 1mm (.040 in) thick, Pre-slit
9-SP1005	9-425	100PCS/PK	9mm Blue PTFE/White Silicone 1mm (.040 in) thick, Pre-slit
9-SP1006	9-425	100PCS/PK	9mm Red PTFE/White Silicone/Red PTFE 1mm (.040") thick
9-SP1007	9-425	100PCS/PK	9mm Red PTFE/White Silicone 1mm (.040 in) thick, Pre-slit (Mercedes Shape)
9-SP1008	9-425	100PCS/PK	9mm White PTFE/Red Silicone 1mm (.040 in) thick, Pre-slit
9-SP1009	9-425	100PCS/PK	9mm Transparent PTFE/Natural Rubber red-orange Septa 1mm (.040 in) thick
9-SP1010	9-425	100PCS/PK	9mm Natural PTFE/White Silicone 1mm (.040") thick, Pre-slit



White PTFE/Red Silicone



Blue PTFE/White Silicone

6mm Micro Insert

P/N	Size	Packing	Description
2ML-N2001	5.8x31.5	100PCS/PK	Insert, Clear Glass, Flat Base
2ML-N2002	5.8x28.5	100PCS/PK	Insert, Clear Glass, Conical Base with Polyspring.
2ML-N2004	5.8x31.5	100PCS/PK	Insert, Clear Glass, Conical Base
SN2004	5.8x28.5	100PCS/PK	PP Insert, Conical Base

- Alllabware Scientific also supply Micro Inserts Vials, including 6mm (5.8mm) inserts. Inserts vials are an economical choice for many routine HPLC applications.
- Polyspring inserts are self-centering with a spring to provide a cushion against needle contact.

