



Chemical Storage Vials

- Storage vials also named sample vials, used in the package of biological agent, cosmetics, high value chemistry and so on.
- Capacity range from 2ml to 60ml.
- Provide constant pH for duration of sample storage life.
- Chemically inert can be suitable for most chromatography and storage applications.



P/N	Size	Packing	Description
2ML-9-V1001C	12x32	100PCS/PK	2ml, Clear Glass, Storage Vial
2ML-9-V1003C	12x32	100PCS/PK	2ml, Amber Glass, Storage Vial
4ML-13-V1001C	15X45	100PCS/PK	4ml, Clear Glass, Storage Vial
4ML-13-V1003C	15X45	100PCS/PK	4ml, Amber Glass, Storage Vial
7ML-15-V1001C	22x40	100PCS/PK	7ml, Clear Glass, Storage Vial
7ML-15-V1003C	22x40	100PCS/PK	7ml, Amber Glass, Storage Vial
8ML-15-V1001C	17x60	100PCS/PK	8ml, Clear Glass, Storage Vial
8ML-15-V1003C	17x60	100PCS/PK	8ml, Amber Glass, Storage Vial
10ML-18-V1001C	22x52	100PCS/PK	10ml, Clear Glass, Storage Vial
10ML-18-V1003C	22x52	100PCS/PK	10ml, Amber Glass, Storage Vial
16ML-18-V1001C	22x70	100PCS/PK	16ml, Clear Glass, Storage Vial
16ML-18-V1003C	22x70	100PCS/PK	16ml, Amber Glass, Storage Vial
20ML-24-V1001C	27.5x57	100PCS/PK	20ml, Clear Glass, Storage Vial
20ML-24-V1003C	27.5x57	100PCS/PK	20ml, Amber Glass, Storage Vial
30ML-24-V1001C	27.5x84	100PCS/PK	30ml, Clear Glass, Storage Vial
30ML-24-V1003C	27.5x84	100PCS/PK	30ml, Amber Glass, Storage Vial
40ML-24-V1001C	27.5x95	100PCS/PK	40ml, Clear Glass, Storage Vial
40ML-24-V1003C	27.5x95	100PCS/PK	40ml, Amber Glass, Storage Vial
60ML-24-V1001C	27.5x140	100PCS/PK	60ml, Clear Glass, Storage Vial
60ML-24-V1003C	27.5x140	100PCS/PK	60ml, Amber Glass, Storage Vial



Cap & Septa

P/N	Size	Packing	Description
9-SP2001-1C	9-425	100PCS/PK	9-425 Blue, Smooth Closed PP Cap
13-SP2001C-B	13-425	100PCS/PK	13-425 Black, Closed Top PP Cap
13-SP2001C-W	13-425	100PCS/PK	13-425 White, Closed Top PP Cap
13-SP2001C-R	13-425	100PCS/PK	13-425 Red, Closed Top PP Cap
13-SP2001C-G	13-425	100PCS/PK	13-425 Green, Closed Top PP Cap
15-SP2001C-B	15-425	100PCS/PK	15-425 Black, Closed Top PP Cap
15-SP2001C-W	15-425	100PCS/PK	15-425 White, Closed Top PP Cap
18-SP2001C-B	18-400	100PCS/PK	18-400 Black, Closed Top PP Cap
18-SP2001C-W	18-400	100PCS/PK	18-400 White, Closed Top PP Cap
20-SP2001C-B	20-400	100PCS/PK	20-400 Black, Closed Top PP Cap
20-SP2001C-W	20-400	100PCS/PK	20-400 White, Closed Top PP Cap
22-SP2001C-B	22-400	100PCS/PK	22-400 Black, Closed Top PP Cap
22-SP2001C-W	22-400	100PCS/PK	22-400 White, Closed Top PP Cap
24-SP2001C-B	24-400	100PCS/PK	24-400 Black, Closed Top PP Cap
24-SP2001C-W	24-400	100PCS/PK	24-400 White, Closed Top PP Cap
28-SP2001C-B	28-400	100PCS/PK	28-400 Black, Closed Top PP Cap
28-SP2001C-W	28-400	100PCS/PK	28-400 White, Closed Top PP Cap
33-SP2001C-B	33-400	100PCS/PK	33-400 Black, Closed Top PP Cap
33-SP2001C-W	33-400	100PCS/PK	33-400 White, Closed Top PP Cap
38-SP2001C-B	38-400	100PCS/PK	38-400 Black, Closed Top PP Cap
38-SP2001C-W	38-400	100PCS/PK	38-400 White, Closed Top PP Cap
45-SP2001C-B	45-400	100PCS/PK	45-400 Black, Closed Top PP Cap
45-SP2001C-W	45-400	100PCS/PK	45-400 White, Closed Top PP Cap
53-SP2001C-B	53-400	100PCS/PK	53-400 Black, Closed Top PP Cap
53-SP2001C-W	53-400	100PCS/PK	53-400 White, Closed Top PP Cap



P/N	Size	Packing	Description
9-SP1001C	9-425	100PCS/PK	9mm Red PTFE/White Silicone 1mm(.040")thick
13-SP1001C	13-425	100PCS/PK	13mm Red PTFE/White Silicone 1.5mm(.060")thick
15-SP1002C	15-425	100PCS/PK	15mm Natural PTFE/White Silicone 1.5mm(.060")thick
18-SP1003C	18-400	100PCS/PK	16.5mm Natural PTFE/White Silicone 1.5mm(.060")thick
20-SP1001C	20-400	100PCS/PK	18.5mm Natural PTFE/White Silicone 1.5mm(.060")thick
22-SP1001C	22-400	100PCS/PK	21.5mm Natural PTFE/White Silicone 1.5mm(.060")thick
24-SP1003C	24-400	100PCS/PK	22mm Natural PTFE/White Silicone 1.5mm(.060")thick
28-SP1001C	28-400	100PCS/PK	27mm Natural PTFE/White Silicone 3.0mm(.120")thick
33-SP1001C	33-400	100PCS/PK	31mm Natural PTFE/White Silicone 3.0mm(.120")thick
38-SP1001C	38-400	100PCS/PK	37.5mm Natural PTFE/White Silicone 1.5mm(.060")thick
45-SP1001C	45-400	100PCS/PK	43mm Natural PTFE/White Silicone 1.5mm(.060")thick
53-SP1001C	45-400	100PCS/PK	53mm Natural PTFE/White Silicone 1.5mm(.060")thick

