

## Screw Top Headspace Vials & Closures

- Screw headspace vials, caps and septa are designed especially for use with Agilent's and other rotating or robotic arm samplers .



### Screw Top Headspace Vial

| P/N           | Size      | Packing   | Description                             |
|---------------|-----------|-----------|---|
| 10ML-18-V1001 | 22.5x46   | 100PCS/PK | 10ml, Clear Glass, Screw Headspace Vial |
| 20ML-18-V1001 | 22.5x75.5 | 100PCS/PK | 20ml, Clear Glass, Screw Headspace Vial |



### Screw Top Headspace Cap & Septa

- Caps are made of magnetic precision metal to exact manufacturing tolerances and lined in a controlled manufacturing environment .
- Pre-slit septa are easier to pierce with needles.
- Pre-assembled caps and septa are convenient and minimize contamination from handling.

| P/N       | Size     | Packing   | Description                               |
|-----------|----------|-----------|---|
| 18-SP2001 | 18mm dia | 100PCS/PK | 18mm Silver Open Top metal Cap (8mm hole) |

| P/N       | Size | Packing   | Description                                |
|-----------|------|-----------|--|
| 18-SP1001 | 17.5 | 100PCS/PK | Blue PTFE/White Silicone 1.5mm(.060")thick |
| 18-SP1002 | 17.5 | 100PCS/PK | Red PTFE/White Silicone 1.5mm(.060")thick  |

