



H1850 Tabletop High Speed centrifuge

10 rotors and many kinds of adaptors are available for H1850, 6x50ml(round/conical) are widely used among them.It has high speed rotor (18500rpm)and high capacity rotor(4x100ml).The special wind deflector can keep the sample in a good state,and the damping system can reduce noise(below 60dB(A)) .



Technical Parameters:

Model	H1850	Noise	≅ 60dB (A)
Order No.	249000	Power Supply	AC220V 50/60Hz 10A
Max.speed	18500rpm	Wattage	450w
Max.RCF	23797xg	Dimension(WxDxH)	420x530x370mm
Max.Capacity	4x100ml	Package Dimension(WxDxH)	487x762x422mm
Speed Accuracy	±50rpm	Net Weight	48kg
Time setting range	1min to 99min59s		

Rotor Options:



Order No.: 805000
 Max.Speed: 18500rpm
 Max.RCF: 23797xg
 RotorCap.: 12x1.5/2.2ml
 Adaptor: 805004(0.5ml)
 805005(0.2ml)



Order No.: 806200
 Max.Speed: 12000rpm
 Max.RCF: 13200xg
 RotorCap.:8x15mlfalcon
 Adaptor: 806205(5ml)



Order No.: 805300
 Max.Speed: 12500rpm
 Max.RCF:16060xg
 RotorCap:12x10/5ml
 Adaptor: 805306(5ml)



Order No.: 805400
 Max.Speed: 11000rpm
 Max.RCF:12850xg
 Rotor Cap.: 6x50ml
 Adaptor: 805405(30ml)
 805406(20ml)
 805408(10ml)
 805409(5ml)
 805410(1.5ml)
 805412(25ml)



Order No.: 805900
 Max.Speed: 11000rpm
 Max.RCF: 12580xg
 RotorCap: 6x50mlfalcon



Order No.: 805500
 Max.Speed: 14000rpm
 Max.RCF:17510xg
 Rotor Cap:24x1.5/2.2ml
 Adaptor: 805005(0.2ml)
 805004(0.5ml)



Order No.: 807000
 Max.Speed: 13000rpm
 Max.RCF: 16410xg
 Rotor Cap.:48x0.5ml



Order No.: 805700
 Max.Speed: 10000rpm
 Max.RCF: 10744xg
 Rotor Cap.:4x100ml



Order No.: 852300
 Max.Speed: 4000rpm
 Max.RCF: 2200xg
 Max.Cap.:2x2x96well



Order No.: 852351
 microplate bucket
 Max.Cap.:2x2x96well



Order No.: 852400
 Max.Speed: 5000rpm
 Max.RCF: 3773xg
 Max.Cap.:4x100ml



Order No.: 852410
 Round bucket
 Cap.:4x100ml

Adaptors:
 852412(50ml falcon)
 852414(50ml)
 852411(15ml falcon)
 852413(15ml)