



HUITAI Drug Stability Test Chamber



◆ Safety Functions

- High and Low Temperature Alarm Emitted by Sound and Light

◆ Recording Temperature and Humidity Data

- Recording the Temperature and Humidity Data through the Built-in printer or USB interface during the stability test running

◆ Features

- Microprocessor Controlling Temperature with colour display
- Temperature Range : 0 - 65 degc
- Humidity Range : 40-95%RH / +/- 3 %
- Can programmable the test running
- Chamber Size provided : 80 L or 150 L or 250 L
- 304 Stainless Steel for Chamber Material
- Two set of Compressors for alternative operation
- With Observable Window
- Touch Screen Display (Optional)

◆ Protections Included :

- Compressor Overheated Protection , Fan Overheated Protection
- Over Temperature Protection, Compressor Over -pressure Protection
- OverLoad protection, Empty water in Water Tank Protection